



## A Comprehensive Solution for HPV Genotyping

Distinguish 28 types of HPV in a single reaction





14 high-risk HPV genotypes

14 low-risk HPV genotypes

**HPV** is the most common sexually transmitted infection (STI). In 2022, approximately 831,204 cancer cases (11.4% in men, 88.6% in women) and 422,935 deaths (10.4% in men, 89.6% in women) were attributed to HPV infection, with cervical cancer accounting for the largest proportion (75.6% of cases).


Age-standardized incidence rate of eight HPV-related cancers by geographic region.



Zhang, J., Ke, Y., Chen, C., et al. HPV cancer burden by anatomical site, country, and region in 2022[J]. Scientific Report, (2025) 15: 21048

## Technology


Patented **Multicolor Melting Curve Analysis (MMCA®)** and **MeltArray** maximize the number of targets per channel, delivering high-level multiplex detection in a single reaction while saving samples and time.




MeltArray 2D Label = Fluorescence Channels + Melting Curve

1 test, 28 HPV genotypes

14 High-Risk HPV



14 Low-Risk HPV



## Advantages



**28 HPV genotypes** identified in a single reaction



**94 samples** in **2.5 hours**



Unique **MMCA®/ MeltArray** multiplex PCR to save precious samples



**Automated** workflow from sample to result

# Specifications

Items	Information
Types	14 High-Risk HPV: 16, 18, 31, 33, 35, 39, 45, 51, 52, 56, 58, 59, 66, 68 14 Low-Risk HPV: 6, 11, 26, 40, 42, 43, 44, 53, 54, 61, 69, 70, 72,
Reagent type	Pre-loaded lyophilized reagent
Sample size	25 µL
Detection specificity	>99%
LoD	500 copies/reaction
Reference gene	GAPDH
Internal control	Positive and negative control are included



## HPV Genotyping Detection Workflow

### Option 1 → Integrated Workflow



#### Sanity 2.0s All-in-One Platform

- Sample-in, result-out
- Minimal hands-on time
- Pre-loaded reagents
- Small footprint

(Flexible loading of 1-4 samples)

### Option 2 → High-Throughput Detection Workflow



**Automated Nucleic Acid Extraction**  
(Up to 96 samples)



**qPCR Detection**  
(Up to 94 samples)




**Data Interpretation**



# WHO guideline for screening and treatment of cervical pre-cancer lesions for cervical cancer prevention

## ● Screening for 14 types of high-risk HPV



## ● Verification of experimental results

was originally found to be 56%, but, when adjusted, this accuracy rose to 80% (Table 2).

### DISCUSSION

Near-patient hrHPV testing in LMICs has the potential to improve outcomes of cervical cancer screening programs by reducing screening costs with use of high-throughput PCR devices, improving care efficiency, and improving retention of at-risk participants. The ideal LMIC test should be sensitive, specific, user friendly, and robust; require minimal equipment; and be rapid and

affordable.<sup>31</sup> With those attributes in mind, hrHPV MMCA near-patient testing was implemented in two separate field locations, each providing representative challenges typical of LMICs.

The hrHPV MMCA assay was highly sensitive and specific at both field sites and accuracy of HPV genotyping enabled additional participant risk stratification, potentially reducing the impact of screening on a taxed health-care system. Since the link of cervical cancer to HPV infection was made, approximately 70% of cervical cancer cases have been linked to infection with HPV16 or

**Table 2.** Genotyping Accuracy of Multicolor Melt Curve Analysis Using Field Protocol

HPV Type	Accuracy (%)	Adjusted Accuracy (%)
16	100	100
18	100	100
59	78	78
68	88	100
58	91	100
31	100	100
39	100	100
Confections	56	80

NOTE: Out-of-the-box multicolor melt curve analysis was used to

box tests to avoid multiple procedural steps. The hrHPV MMCA achieves this by providing lyophilized hrHPV master mix prepackaged into PCR strip tubes. Temperatures at field sites in Honduras reached > 32.2°C, with high levels of humidity and no access to refrigeration for reagent storage. After 4 days of exposure to elevated temperatures, no differences in test sensitivity or specificity were observed. In addition, minimally trained personnel needed only to pipet crude lysate directly into the lyophilized reagents and place the reaction tubes into the quantitative PCR (qPCR) thermocycler for successful analysis. As noted, invalid rates were highest during the testing at the first field site but dropped in the more infrastructure-challenged

**Zeesan HPV Assay** offers high sensitivity and specificity, which reduces the risk of progression to cervical precancer / cancer due to undefined genotypes.

Experimental results shows the high stability of **Zeesan HPV Assay**.

Turner S A , Deharvenst S J , Lyons K D , et al. Implementation of Multicolor Melt Curve Analysis for High-Risk Human Papilloma Virus Detection in Low- and Middle-Income Countries: A Pilot Study for Expanded Cervical Cancer Screening in Honduras[J]. Journal of Global Oncology, 2018(4):1-8.

## Product Information

Cat.No.	Product name	Size	Compatible Platform
608112	HPV DNA Extraction Kit (Sanity 2.0)	24 Tests/Kit	Sanity 2.0s System
801170	High Risk HPV Genotyping Kit (Sanity 2.0)	24 Tests/Kit	
801946	MeltArray High Risk HPV Genotyping Assay (Sanity)	24 Tests/Kit	
801916	Low Risk HPV Genotyping Kit (Sanity 2.0)	24 Tests/Kit	
604106	Lab-Aid 824 HPV DNA Extraction Kit	48 Tests/Kit	824s/896 Extraction System
610104	HPV DNA Extraction Kit (Lab-Aid 896)	96 Tests/Kit	
801011	MeltPro® High Risk HPV Genotyping Assay (MMCA)	48 Tests/Kit	
801943	MeltArray High Risk HPV Genotyping Assay	48 Tests/Kit	SLAN Real-time PCR System
801013	MeltPro® Low Risk HPV Genotyping Assay	48 Tests/Kit	

For more information, please visit <http://www.zeesandx.com>

**Xiamen Zeesan Biotech Co., Ltd.**

No. 884-1 and 884-2 Lianting Road,  
(Xiang An) Industrial Area, Torch High-Tech Zone,  
Xiamen City, Fujian Province, 361101, P.R. China

**Phone:** +86-592-7615091

**E-mail:** info@zsandx.com

**Website:** www.zeesandx.com

